INDEX, ÍNDICE, INDICE, ÍNDICE	
2 – 7	English
8 – 14	Español
15 – 20	Italiano
21 – 27	Français
28 – 33	Português
34 – 35	Symbols Símbolos Simboli Symboles Símbolos

English

VersaTREK® MYCO PZA KIT, 7116-70

INTENDED USE

The VersaTREK Myco PZA Kit is a rapid qualitative procedure for *in vitro* susceptibility testing of *Mycobacterium tuberculosis*, from culture, to pyrazinamide (PZA). The VersaTREK Myco PZA Kit is used with the VersaTREK Automated Microbial Detection System (VTI) and the ESP® Culture System II (ESP).

Appropriate inoculum sources are prepared from a liquid VersaTREK Myco seed bottle.

SUMMARY AND EXPLANATION

Rapid diagnosis of tuberculosis is important for initiation of appropriate drug therapy and institution of measures to prevent the spread of this highly contagious disease. Delayed treatment, particularly of multiple-drug-resistant *M. tuberculosis* (MDR-TB) strains, can result in treatment failure and patient mortality.¹ To address these issues, clinical microbiologists are urged to use the most sensitive and rapid diagnostic tests that provide reliable results. Experts at the CDC have made several recommendations regarding Mycobacterial testing, including the use of both a liquid and solid medium for Mycobacterial culture.²

For pyrazinamide susceptibility testing of *M. tuberculosis*, the VTI/ESP System combines a liquid culture medium (VersaTREK Myco Broth), a growth supplement (VersaTREK Myco GS), and a specific concentration of pyrazinamide with a detection system that automatically incubates and continuously monitors culture bottles inoculated with isolates of *M. tuberculosis* from a variety of specimen sources. The sponges in the VersaTREK Myco bottles provide a growth support matrix and increase the surface area exposed to headspace oxygen. The technology of the VTI/ESP System is based on the detection of headspace pressure changes within a sealed bottle due to microbial growth. A special algorithm has been developed for detection of very slow growing Mycobacteria.

With the VTI/ESP System, anti-tuberculosis drug PZA can be injected into VersaTREK Myco bottles supplemented with VersaTREK Myco GS. An inoculum of a confirmed culture of *M. tuberculosis* grown from a seed bottle is added to a drug-free control bottle, used as a positive control, and to a test bottle containing PZA. If the organism is susceptible to the drug, delayed or no microbial growth will be detected. If the organism is resistant to PZA, microbial growth will be detected by the generation of a curve exhibiting a downward excursion from the baseline. The time to detection of drug-free bottle determines when a test is complete. The drug used in this system is the primary drug, pyrazinamide, which is hydrated with an acidic rehydration buffer to a specific critical concentration. The purpose of the acidic rehydration buffer is to attain a final pH of 5.9 to 6.0, which optimizes the activity of the PZA drug. The specific drug concentration for this test is 300 µg/ml.

The detection of microbial growth in the VTI/ESP System is based on the sensitive measurement of gas production and/or gas consumption within the headspace of a sealed culture bottle. This information is used to generate a pressure curve for each bottle. An internal algorithm analyzes the information to determine the status of each test sample. When a set of conditions is met, a bottle is flagged as positive. There are approximately 10⁶ CFU/ml of Mycobacteria at the time of detection.

PRINCIPLES OF THE PROCEDURE

Susceptibility testing of *M. tuberculosis* specimens is performed with the VTI/ESP System by inoculating a suspension of organisms into VersaTREK Myco broth bottles that have been supplemented with VersaTREK Myco GS and pyrazinamide. A VersaTREK Connector is attached to each bottle to establish a closed system for monitoring the bottle with the VTI/ESP instrument. A hydrophobic membrane in the VersaTREK Connector prevents aerosolization.

Isolate information is entered into the system computer and the bottles are placed in the instrument for incubation at 35 – 37 °C under stationary conditions. VTI/ESP System detects Mycobacterial growth by continuously monitoring (every 24 minutes) the rate of consumption of oxygen within the headspace of the culture bottle and reports that growth response with a positive signal. At the end of the specified incubation period, as determined by the drug-free control bottle for each isolate that is tested, the isolate is manually determined to be susceptible or resistant to the drug.

REAGENTS

The VersaTREK Myco PZA Kit contains two bottles of the lyophilized pyrazinamide drug and four bottles of the acidic PZA Rehydration Buffer. Each drug bottle is to be rehydrated to 25 ml with the acidic rehydration buffer.

FORMULA UPON REHYDRATION*

(PRB)	Acidic PZA Rehydration Buffer	25 ml
(PZA)	Pyrazinamide	0.45% w/v
(SA)	Solubilizing Agent	1.25% w/v

*Adjusted to meet performance criteria

PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Refer to Material Safety Data Sheets for complete hazard information.
3. **WARNING! POTENTIAL INFECTIOUS TEST SPECIMEN.** Follow Universal Precautions and institutional policy in handling and disposing of infectious agents.
4. When working with *Mycobacterium tuberculosis* grown in culture, the recommendation is for Biosafety Level 3 practices, containment equipment and facilities.³
5. The VersaTREK Connector contains a sharp recessed needle that is sterile upon first removing the seal, but is potentially contaminated upon termination of use.

Follow institutional policy for handling and disposing of biohazard-contaminated devices. (Do not handle in a casual manner. Sterilize prior to disposal. Dispose in an approved sharps container. Do not reuse.)
6. Visually inspect all bottles for contamination, cracks or other signs of deterioration. Do not use bottles that appear turbid or cracked.
7. Use only needle locking syringes.
8. Autoclave all VersaTREK Myco bottles prior to disposal.
9. Rehydrate the lyophilized PZA drug bottles only with the acidic PZA Rehydration Buffer provided in the kit.

STORAGE

VersaTREK Myco PZA: Store at 2-8°C. Protect from light. Upon rehydration, store at 2-8°C and use within 4 weeks, or freeze at -20°C and use within 3 months. Once thawed, do not refreeze.

PROCEDURE: SUSCEPTIBILITY TESTING OF *MYCOBACTERIUM TUBERCULOSIS* USING THE VersaTREK/ESP CULTURE SYSTEM II

Materials Provided:

VersaTREK Myco Culture Bottles	7111-42
VersaTREK Myco GS	7112-42
VersaTREK Myco PZA Kit	7116-70
VersaTREK Connectors	7150-44

Materials Required but Not Provided:

Sterile distilled or deionized water
Autoclave
Middlebrook 7H10 or other mycobacterial agar or egg-base medium, for purity checks
CO₂ incubator, 37°C
Sterile loops
Sterile tuberculin syringes with permanently attached needles
Mycobactericidal disinfectant (e.g. Amphyl™, OMNI™, etc.)
Alcohol swabs
Vortex mixer
Biological safety cabinet
Microscope
Materials for staining slides
Quality Control organism: *M. tuberculosis* ATCC® 27294

Inoculation of VersaTREK Myco Culture Bottles for Myco Susceptibility Testing

1. Label each VersaTREK Myco culture bottle with identifying information. For PZA testing, 2 bottles will be used per test.
2. Disinfect the culture bottle stoppers with alcohol.
3. Aseptically add 1.0 ml of VersaTREK Myco GS to each culture bottle by injecting through the septum with a needle and syringe.
4. Rehydrate lyophilized PZA by aseptically adding 25 ml acidic PZA Rehydration Buffer to each drug vial. This yields a drug concentration of 4500 µg/ml.

- Aseptically add 1.0 ml of the drug to each culture bottle by injecting through the septum with a needle and syringe. The final concentration of pyrazinamide is 300 µg/ml. For the drug-free control bottle, add 1 ml of the acidic PZA Rehydration Buffer in place of the drug. Twenty-five tests can be performed from a single vial.

Drug	Concentration of Drug after Reconstitution	Volume added to Myco test bottles	Final Concentration in Test bottles
PZA	4500 µg/ml	1 ml	300µg/ml

- The inoculum of a specimen for susceptibility testing must be confirmed as a pure culture of *M. tuberculosis* prior to PZA testing and prepared as follows:

VersaTREK Myco Seed Bottle - A 0.5 ml volume of a 1.0 McFarland equivalent created from an agar or liquid source (growth in Middlebrook 7H9 or liquid from a primary detection bottle) is aseptically injected into a VersaTREK Myco culture bottle supplemented with 1.0 ml VersaTREK Myco GS. A sterile VersaTREK Connector is placed on the bottle and the bottle is incubated in the instrument until a positive signal is determined. This seed bottle must be used within 72 hours. Remove the VersaTREK Connector and vortex for 1-2 minutes. Add 0.5 ml of a 1:10 dilution of culture from a positive VersaTREK Myco seed bottle to the appropriate culture bottles containing drug and to the drug-free growth control bottle.

- Wipe the culture bottles and stoppers with a Mycobactericidal disinfectant.
- Mix all bottles by inverting several times.
- Remove the seal from the bottom of a VersaTREK Connector. Press the needle end of the VersaTREK Connector through the culture bottle stopper. NOTE: Do not invert the culture bottle while the connector is on the bottle. Fluid in the needle may interfere with the pressure readings for the bottle.
- Place the culture bottle with the VersaTREK Connector into the instrument. Refer to the VTI/ESP System Operator's Manual for details on bottle entry.

Loading suggestion: If the sample accession number is "123", give the PZA test bottles accession numbers such as 123PC for the positive growth control, and 123P for pyrazinamide. In this way, all bottles associated with this sample accession number will be able to be searched for and located as a group.

RESULTS

For susceptibility testing using the VTI/ESP System, a test isolate is interpreted as being susceptible or resistant to a drug based on a formula. Record the time to detection of the drug-free bottle. This time is used to calculate susceptible or resistant interpretations using the following formula:

Susceptible: No growth, or the time to detection of the drug bottle is greater than three days (rounded to the nearest whole number) of the time to detection of the drug-free control bottle (rounded to nearest whole number). To round to the nearest whole number, if time ≤ XX.49 is round down and if time ≥ XX.50 is round up.

Resistant: Time to detection of drug bottle is ± three days (rounded to the nearest whole number) of the time to detection of the drug-free control bottle (rounded to the nearest whole number). To round to the nearest whole number, if time ≤ XX.49 is round down and if time ≥ XX.50 is round up.

NOTE:

- No bottle should flag positive before the third day. Bottles that flag before three days are either contaminated or the inoculum is too high. The test should be repeated.
- A test bottle may signal before the control bottle (see Table 2 Example D).

All resistant results should be confirmed by an alternate method prior to reporting.

Bottles that signal positive should be stained for AFB (Acid Fast Bacilli) and subcultured to a Middlebrook 7H10 or 7H11 agar plate, as well as a non-selective agar plate (e.g. Blood Agar) to ensure purity of the specimen. Growth control bottles that signal in less than 2.5 days may indicate the presence of a non-Mycobacterial contaminant whereas control bottles that do not signal within 12 days after being inoculated may be indicative of improper bottle preparation. For either situation, the specimen should be retested.

Mono-resistance to PZA is uncommon, therefore in the event of unexpected resistance results, verify purity and identification of the isolate.

ALWAYS review the isolate's graph before reporting. Refer to the Operator's Manual for graph analysis. If graphs for tests that have signaled positive look abnormal, contact TREK Technical Support at 1-800-642-7029 or 1-216-351-8735.

PZA Result reporting:

TABLE 1: Test Observations and Test Assessment

Test bottle	Result obtained	Possible cause	Action
Growth control	<3 days from inoculation	Possible contaminant or too heavy inoculum; check purity plates and graphs.	REPEAT TEST
Growth control*	>12 days from inoculation	Non viable organism, too light inoculum	REPEAT TEST
Growth control	Growth detected between 3 and 12 days after inoculation		Valid test - evaluate drug bottles
Purity check	No contaminants observed by at least 18 h on Blood Agar plate		Valid test
	Non-AFB organisms observed	Culture was not pure	REPEAT TEST
	No colonial variation observed on Middlebrook agar plates		Valid test- No action
	Colonial variation observed on Middlebrook agar plates	More than one species of mycobacterium is present	REPEAT TEST to verify results; notify physician that final results are pending

*The upper range of 12 days detection time for the growth control was established by clinical and in-house data.

Be sure to check the subculture results on the mycobacterial agar to verify that only one species of mycobacterium was present, especially in the case of a resistant result.

TABLE 2: Examples of Results

Patient	A	B	C	D	E	F	G
PZA Control Positive (Days)	2	4	5	9	10	12	13
PZA Drug Positive (Days)	5	2	6	7	14	No Growth	14
Result	Invalid Growth Control, Repeat Test	Invalid, Repeat Test	Resistant*	Resistant*	Susceptible	Susceptible	Invalid Growth Control, Repeat Test

*All resistant results should be confirmed by an alternate method prior to reporting.

QUALITY CONTROL

A new lot number or shipment of VersaTREK Myco PZA should be tested with the quality control organism *M. tuberculosis* ATCC 27294 (PZA susceptible). Expected results are shown in Table 3. Quality control should be run once each week when susceptibility testing is performed. Observation of the proper results, within 4-15 days, indicates the reagents are ready for use in testing patient isolates.

If proper results are not observed, do not report patient's results. Repeat quality control and any patient isolates affected by the initial quality control failure. If repeat quality control does not perform as expected, do not report patient results. Do not use the product until you have contacted Technical Services.

1. Add 1 ml VersaTREK Myco GS to VersaTREK Myco culture bottles.
2. Rehydrate lyophilized pyrazinamide with acidic PZA Rehydration Buffer as described under "Inoculation of VersaTREK Myco Culture Bottles for Susceptibility Testing (Item 4)."
3. Add 1 ml of drug to each VersaTREK Myco bottle to yield final drug concentrations of pyrazinamide at 300 µg/ml. Add 1.0 ml PZA Rehydration Buffer to the growth control bottle.

4. Add 0.5 ml of a 1:10 dilution of culture from a positive VersaTREK Myco seed bottle to the appropriate culture bottles containing drug and to the drug-free growth control bottle.
5. After observation of expected results, use the remainder of the lot in testing specimens. If expected results are not achieved, contact VersaTREK Technical Support.

TABLE 3 Expected Quality Control Results

Organism	ATCC	Pyrazinamide
<i>M. tuberculosis</i>	27294	Susceptible

The most common causes of quality control failure are over-inoculation of the PZA test sets and contamination. If repeated false resistance results are obtained, the PZA antibiotic most likely has lost activity. The recommendation is to test with a new antibiotic stock.

LIMITATIONS

1. The VersaTREK Myco PZA susceptibility test does not interpret the degree of susceptibility of the isolate being tested. It is a qualitative test.
2. The VersaTREK Myco PZA susceptibility test can only be performed on the VTI or ESP instruments. Testing cannot be performed manually.
3. Use only pure cultures of *M. tuberculosis*. The user must confirm the purity of the culture, especially if the result is resistant. Cultures that contain multiple species of Mycobacteria may give erroneous results.
4. Failure to prepare the isolates per the manufacturer's recommendation may lead to inaccurate results.
5. Failure to re-hydrate the PZA with the proper volume of acidic PZA Rehydration Buffer may give inaccurate results.
6. Thorough mixing of inoculated bottles is important. Failure to mix adequately may lead to false resistant results.
7. Failure to use the VersaTREK Myco GS growth supplement may lead to false susceptible results.

EXPECTED RESULTS

In the initial testing, a total of 96 clinical isolates of *M. tuberculosis* were tested with the VersaTREK Myco PZA Kit at five geographically diverse sites. The testing included both fresh clinical and sub-cultured isolates from liquid and solid media sources.

The repeat rates for this study was <10%. The average time to result for the VersaTREK Myco PZA susceptibility test was 8.18 days for a liquid seed inoculum, with a range from 3 to 15 days. See figure 1 below.

PERFORMANCE CHARACTERISTICS

In house testing

Lot Reproducibility:

Lot reproducibility was performed using 4 well-characterized strains of *M. tuberculosis*, which were tested in triplicate and on three separate days. Three lots of VersaTREK Myco broth were represented in the testing. The overall reproducibility was 100% for seeded inoculum.

Eight lots of manufactured VersaTREK Myco PZA, seven VersaTREK Myco broth, and six VersaTREK Myco GS were evaluated for reproducibility with acceptable results for all lots.

CDC Challenge Panel Testing:

Ten well-characterized strains of *M. tuberculosis* were tested and compared to BACTEC 460TB and to expected results. Eight strains were susceptible and two were resistant to PZA. There was an overall 100% agreement with BACTEC and expected results.

Comparison of ESP Culture System II and VersaTREK Microbial Detection System:

Recovery of different microorganisms and time to detection was used to demonstrate no difference in performance between ESP Culture System II and VersaTREK Microbial Detection System with additional testing of *M. tuberculosis* and PZA. Results and time to detection by either system were not significantly different.

Clinical Site testing:

Clinical Evaluation:

The VersaTREK Myco PZA susceptibility test was evaluated at five geographically diverse sites composed of regional reference centers, public health laboratories, and university hospital-based laboratories. The VersaTREK Myco PZA Kit was compared to the BACTEC 460TB PZA susceptibility test method.

Reproducibility Testing:

Lot reproducibility was performed using two well-characterized strains of *M. tuberculosis* from CDC using two lots each of VersaTREK Myco broth and VersaTREK Myco GS. Testing was performed in triplicate and on three separate days. The overall reproducibility was $\geq 95\%$ for seeded inoculum.

Challenge Testing:

CDC challenge set of 30 organisms was tested at three sites. Of the possible 90 test points, only 77 were valid. The lower number of tests (7.2% repeat rate) was due to no growth, contamination or instrument issues. Of the 77 organisms tested, 68 were susceptible and 9 were resistant. There was an overall agreement of 98.7% with expected results

Clinical Isolate Testing using Seed Inocula

Susceptibility tests of *Mycobacterium tuberculosis* were performed on 96 isolates from a liquid inoculum source.


Table 4 presents the results from clinical isolate testing for PZA at 300 $\mu\text{g/ml}$ for ESP and 100 $\mu\text{g/ml}$ for BACTEC. Testing showed 11 discordant results for seeded inoculum. Of the 11 discordant results from seeded inoculum, 10 were resistant for ESP/VTI and susceptible for BACTEC; and one was susceptible for ESP/VTI and resistant for BACTEC.

Table 4. Clinical Isolate Results- VersaTREK Myco PZA susceptibility test compared to BACTEC 460TB susceptibility.

		BACTEC 460TB System		ESP/VTI System			
		Expected Results		Susceptible		Resistant	
Source	# Tests	S	R	# Agree	Category Agreement % (95% CI) ^a	# Agree	Category Agreement % (95% CI)
Liquid	96	83	13	73	88 (79-94.1)	12	92.3 (64-99.8)

a. 95% Confidence Intervals (CI). Lee, Eric W. and Dubin, N. Estimation and Sample Size Considerations for Clustered Binary Responses. 1994. Statistics in Medicine. Vol. 13. 1241-1252.

Figure 1. Total test time for seeded inoculum preparation for VersaTREK Myco PZA Susceptibility testing of *Mycobacterium tuberculosis*

HAZARD INFORMATION

Pyrazinamide

Xn	Harmful
R22	Harmful if swallowed.
S20	When using do not eat or drink.
S36	Wear suitable protective clothing.
S46	If swallowed, seek medical advice immediately and show this container or label.

REFERENCES

1. Stratton, C.W. 1992. The resurgence of tuberculosis. Inf. Dis. Newsl. 11:73-80.
2. Tenover, F.C., J.T. Crawford, R.E. Huebner, L.J. Geiter, C.R. Horsberg, Jr., and R.C. Good. 1993. The resurgence of tuberculosis: Is your laboratory ready? J. Clin. Microbiol. 31:767-770.
3. CDC/National Institutes of Health, Biosafety in Microbiological and Biomedical Laboratories. 4th Edition Atlanta: US Department of Health and Human Services, Public Health Service, CDC and NIH, 1999; DHHS publication no. (CDC) 93-8395.

ATCC is a trademark of the American Type Culture Collection

BACTEC is a trademark of Becton Dickinson and Company

ESP and VersaTREK are trademarks of TREK Diagnostic Systems, Inc.

For additional information call TREK Technical Support at 1.800.642.7029 or 1-216-351-8735.

TREK Diagnostic Systems
982 Keynote Circle, Suite 6
Cleveland, Ohio 44131
USA
800.871.8909
www.trekds.com

TREK Diagnostic Systems
Units 17 – 19, Birches Industrial Estate
East Grinstead
West Sussex, RH19 1XZ
UK

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