<b>INDEX, ÍNDICE, INDICE, ÍNDICE</b>	
<b>2 – 7</b>	<b>English</b>
<b>8 – 13</b>	<b>Español</b>
<b>14 – 19</b>	<b>Italiano</b>
<b>20 – 26</b>	<b>Français</b>
<b>27 – 32</b>	<b>Português</b>
<b>33 – 34</b>	<b>Symbols Símbolos Simboli Symboles Símbolos</b>

## English

### **VersaTREK® Myco Streptomycin, 7120-30**

#### **INTENDED USE**

The VersaTREK Myco Streptomycin test is a rapid qualitative procedure for *in vitro* susceptibility testing of *Mycobacterium tuberculosis*, from culture, to streptomycin. VersaTREK Myco Streptomycin is used with the VersaTREK Automated Microbial Detection System (VTI) and the ESP® Culture System II (ESP).

#### **SUMMARY AND EXPLANATION**

Rapid diagnosis of tuberculosis is important for initiation of appropriate drug therapy and institution of measures to prevent the spread of this highly contagious disease. Delayed treatment, particularly of multiple-drug-resistant *M. tuberculosis* (MDR-TB) strains, can result in treatment failure and patient mortality.<sup>1</sup> To address these issues, clinical microbiologists are urged to use the most sensitive and rapid diagnostic tests that provide reliable results. Experts at the CDC have made several recommendations regarding Mycobacterial testing, including the use of both a liquid and solid medium for Mycobacterial culture.<sup>2</sup>

For streptomycin susceptibility testing of *M. tuberculosis*, The VTI/ESP System combines a liquid culture medium (VersaTREK Myco Broth), a growth supplement (VersaTREK Myco GS), and specific concentrations of streptomycin with a detection system that automatically incubates and continuously monitors culture inoculated with isolates of *M. tuberculosis* from a variety of specimen sources. The sponges in the VersaTREK Myco bottles provide a growth support matrix and increase the surface area exposed to headspace oxygen. The technology of the VTI/ESP System is based on the detection of headspace pressure changes within a sealed bottle due to microbial growth. A special algorithm has been developed for detection of very slow growing Mycobacteria.

With the VTI/ESP System, anti-tuberculosis drug streptomycin can be injected into VersaTREK Myco bottles supplemented with VersaTREK Myco GS. An inoculum of a confirmed culture of *M. tuberculosis* grown from a seed bottle is added to a drug-free control bottle, used as a positive control, and to a test bottle containing streptomycin. If the organism is susceptible to the drug, delayed or no microbial growth will be detected. If the organism is resistant to streptomycin, microbial growth will be detected by the generation of a curve exhibiting a downward excursion from the baseline. The time to detection of drug-free bottle determines when a test is complete. The drug used in this system is streptomycin which is hydrated with sterile water to a specific critical concentration. The specific drug concentrations for this drug are 2 µg/ml and 8 µg/ml.

The detection of microbial growth in the VTI/ESP System is based on the sensitive measurement of gas production and/or gas consumption within the headspace of a seeded culture bottle. This information is used to generate a pressure curve for each bottle. An internal algorithm analyzes the information to determine the status of each test sample. When a set of conditions is met, a bottle is flagged as positive. There are approximately 10<sup>6</sup> CFU/ml of Mycobacteria at the time of detection.

#### **PRINCIPLES OF THE PROCEDURE**

Susceptibility testing of *M. tuberculosis* specimens is performed with the VTI/ESP System by inoculating a suspension of organisms into VersaTREK Myco broth bottles that have been supplemented with VersaTREK Myco GS and VersaTREK Myco Streptomycin. A VersaTREK Connector is attached to each bottle to establish a closed system for monitoring the bottle with the VTI/ESP instrument. A hydrophobic membrane in the VersaTREK Connector prevents aerosolization.

Isolate information is entered into the system computer and the bottles are placed in the instrument for incubation at 35 – 37 °C under stationary conditions. VTI/ESP System detects Mycobacterial growth by continuously monitoring (every 24 minutes) the rate of consumption of oxygen within the headspace of the culture bottle and reports that growth response with a positive signal. At the end of the specified incubation period, as determined by the drug-free control bottle for each isolate that is tested, the isolate is manually determined to be susceptible or resistant to the drug.

#### **REAGENTS**

VersaTREK Myco Streptomycin contains one vial of lyophilized streptomycin. The vial is to be rehydrated with 25 ml of sterile distilled or deionized water.

#### **FORMULA UPON REHYDRATION \***

(H2O)	Processed Water	25 ml
(STR)	Streptomycin	0.012% w/v
(SA)	Solubilizing Agent	1.25% w/v

\*Adjusted to meet performance criteria.

## PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Refer to Material Safety Data Sheets for complete hazard information.
3. **WARNING! POTENTIAL INFECTIOUS TEST SPECIMEN.** Follow Universal Precautions and institutional policy in handling and disposing of infectious agents.
4. When working with *M. tuberculosis* grown in culture, the recommendation is for Biosafety Level 3 practices, containment equipment and facilities<sup>3</sup>
5. The VersaTREK Connector contains a sharp recessed needle that is sterile upon first removing the seal, but is potentially contaminated upon termination of use.
6. Visually inspect all bottles for contamination, cracks or other signs of deterioration. Do not use bottles that appear turbid or cracked.
7. Use only needle locking syringes.
8. Autoclave all VersaTREK Myco bottles prior to disposal.

## STORAGE

Store at 2-8°C. Protect from light. Upon rehydration, store at 2-8°C and use within 4 weeks, or freeze at -20°C and use within 3 months. Once thawed, do not refreeze.

## EXPIRATION DATE

The expiration date applies to the product in its intact (unopened) container when stored as directed. The VersaTREK Streptomycin Test should not be used beyond the indicated expiration date.

## PROCEDURE: SUSCEPTIBILITY TESTING OF *M. tuberculosis* USING THE VERSATREK/ESP CULTURE SYSTEM II

### Materials Provided:

VersaTREK Myco Streptomycin	7120-30
VersaTREK Myco Culture Bottles	7111-42
VersaTREK Myco GS	7112-42
VersaTREK Connectors	7150-44

### Materials Required but Not Provided:

Sterile distilled or deionized Water  
Autoclave  
Middlebrook 7H10 or other mycobacterial agar or egg-base medium, for purity checks  
CO<sub>2</sub> incubator, 37°C  
Sterile loops  
Sterile tuberculin syringes with permanently attached needles  
Mycobacterial disinfectant (e.g. Amphly™, OMNI™, etc.)  
Alcohol swabs  
Vortex mixer  
Biological safety cabinet  
Microscope  
Materials for staining slides  
Quality Control Organism; *M. tuberculosis* ATCC® 25618

### Inoculation of VersaTREK Myco Culture Bottles for Myco Susceptibility Testing

1. Label each VersaTREK Myco culture bottle with identifying information. For streptomycin testing, 3 bottles will be used per test.
2. Disinfect the culture bottle stoppers with alcohol.
3. Aseptically add **1.0 ml** of VersaTREK Myco GS to each culture bottle by injecting through the septum with a needle and syringe.
4. Rehydrate lyophilized VersaTREK Myco Streptomycin by aseptically adding 25 ml of sterile distilled or deionized water. This yields a drug concentration of 120 µg/ml.
5. Aseptically remove 5 ml of the solution and transfer into a sterile container. Add 15 ml of distilled or deionized sterile water. This yields a drug concentration of 30 µg/ml.

6. Aseptically add **1.0 ml** of the 120 µg/ml drug concentration to a culture bottle by injecting through the septum with a needle and syringe. The final streptomycin concentration in the Myco culture bottle is 8 µg/ml.
7. Aseptically add **1.0 ml** of the 30 µg/ml drug concentration to a second culture bottle by injecting through the septum with a needle and syringe. The final streptomycin concentration in the Myco culture bottle is 2 µg/ml.

**TABLE 1: Anti-Tuberculosis Drug Preparation**

Drug and Final Concentration in VersaTREK Myco Bottles	Amount of Diluent Added to Lyophilized Streptomycin	Concentration After Rehydration	Additional Dilutions	Concentration to be Injected into the Test Bottles	Number of Tests per Vial
Streptomycin 8.0 µg/ml	25 ml	120 µg/ml	No	120 µg/ml	20
Streptomycin 2.0 µg/ml		↕⇒	5 ml diluted drug plus 15 ml water	30 µg/ml	20

8. For the drug-free control bottle, add **1.0 ml** of sterile water in place of the drug. Twenty tests can be performed from a single vial.
9. The inoculum of a specimen for susceptibility testing must be confirmed as a pure culture of *M. tuberculosis* prior to streptomycin testing and prepared as follows:
  - a) **McFarland equivalent cell suspension** - Pure isolates (not more than 4-5 weeks old) growing on an appropriate mycobacterial agar medium (Middlebrook 7H10, Middlebrook 7H11, Lowenstein-Jensen slants) are scraped using a sterile loop and placed in a tube containing approximately 9 to 10 ml sterile distilled or deionized water and glass beads. The tube is vortexed and allowed to settle for 15 minutes. Approximately 5 ml of the top layer is removed, placed in a sterile test tube and adjusted to a No. 1 McFarland equivalent (approximately  $3.0 \times 10^7$  CFU/ml) using sterile distilled or deionized water.
  - b) **VersaTREK Myco seed bottle** - Prepare a No. 1 McFarland equivalent from an agar or liquid source (growth in Middlebrook 7H9) of the test isolate. Aseptically inject 0.5 ml of this dilution into a VersaTREK Myco culture bottle supplemented with 1.0 ml VersaTREK Myco GS. A sterile VersaTREK Connector is placed on the bottle and the bottle is incubated in the instrument until a positive signal is observed (approx.  $10^6$  CFU/ml). Once a positive signal is reached, the VersaTREK Connector is removed and the bottle is vortexed for 1-2 minutes. The seed bottle inoculum must be used within 72 hours of the positive instrument signal.
10. Perform a 1:10 dilution of the McFarland No. 1 from agar or the Myco seed bottle in sterile saline or distilled water (inoculum for test). Add 0.5 ml to each control and test broth.
11. Wipe the culture bottles and stoppers with a Mycobactericidal disinfectant.
12. Mix all bottles by inverting several times.
13. Remove the seal from the bottom of a VersaTREK Connector. Press the needle end of the VersaTREK Connector through the culture bottle stopper. NOTE: Do not invert the culture bottle while the connector is on the bottle. Fluid in the needle may interfere with the pressure readings for the bottle.
14. Place the culture bottle with the VersaTREK Connector into the instrument. Refer to the VersaTREK/ESP System Operator's Manual for details on bottle entry.
 

Loading suggestion: If the sample accession number is "123", give the streptomycin test bottles accession numbers such as 123PC for the positive control, and 123S2 and 123S8 for the streptomycin test broths. In this way, all bottles associated with this sample accession number will be able to be searched for and located as a group.
15. An alternative method to adding reagents and inocula to broth bottles is to remove the cap and stopper and perform the additions using a sterile pipette. Careful closure of the stopper and cap is required.

## RESULTS

For susceptibility testing using the VersaTREK/ESP Culture System II, an isolate is interpreted as being susceptible or resistant to a drug based on the following (round all times of detection to the nearest whole day):

1. After the drug-free bottle has become positive, record the time to detection (this is found on the graph of the test bottle).
2. Three additional days, rounded to the nearest whole number, of incubation after the drug-free control has signaled positive, search and record the results from the other test bottles for the isolate.
3. **Susceptible:** An isolate in a drug-containing bottle that has not signaled within the three-day detection period of the drug-free control is determined to be susceptible at that drug concentration. Growth after the three days, rounded to the nearest whole number, of the growth control bottle signaling positive should be disregarded because it is not considered to be an indication of resistance. To round to the nearest whole number, if time  $\leq$  XX.49 is round down and if time  $\geq$  XX.50 is round up.
4. **Resistant:** Time to detection of drug bottle is  $\pm$  three days (rounded to the nearest whole number) of the time to detection of the drug-free control bottle (rounded to the nearest whole number). To round to the nearest whole number, if time  $\leq$  XX.49 is round down and if time  $\geq$  XX.50 is round up.

### NOTE:

- No bottle should flag positive before the third day. Bottles that flag before three days are either contaminated or the inoculum is too high. The test should be repeated.
- A test bottle may signal before the control bottle (see Table 3 example D).

**All resistant results should be confirmed by an alternate method prior to reporting.**

Bottles that signal positive should be stained for AFB (Acid Fast Bacilli) and subcultured to a Middlebrook 7H10 or 7H11 agar plate, as well as a non-selective agar plate (e.g. Blood Agar) to ensure purity of the specimen. Growth control bottles that signal in less than 2.5 days may indicate the presence of a non-Mycobacterial contaminant whereas control bottles that do not signal within 10 days after being inoculated may be indicative of improper bottle preparation. For either situation, the specimen should be retested.

Always review the isolate's graph before reporting. Refer to the Operator's Manual for graph analysis. If graphs for tests that have signaled positive look abnormal, contact TREK Technical Support at 1-800-642-7029 or 1-216-351-8735.

**TABLE 2: Test Observations and Test Assessment**

Test bottle	Result obtained	Possible cause	Action
Growth control	<3 days from inoculation	Possible contaminant or too heavy inoculum; check purity plates and graphs.	REPEAT TEST
Growth control*	>10 days from inoculation	Non viable organism, too light inoculum	REPEAT TEST
Growth control	Growth detected between 3 and 10 days after inoculation		Valid test - evaluate drug bottles
Purity check	No contaminants observed by at least 18 h on Blood Agar plate		Valid test
	Non-AFB organisms observed	Culture was not pure	REPEAT TEST
	No colonial variation observed on Middlebrook agar plates		Valid test- No action
	Colonial variation observed on Middlebrook agar plates	More than one species of mycobacterium is present	REPEAT TEST to verify results; notify physician that final results are pending

\*The upper range of 10 days detection time for the growth control was established by clinical and in-house data.

Be sure to check the subculture results on the mycobacterial agar to verify that only one species of mycobacterium was present, especially in the case of a resistant result.

**TABLE 3: Examples of Results for Calculation of Susceptibility of *M. tuberculosis* strains to VersaTREK Myco Streptomycin**

Patient	A	B	C	D	E	F	G
Streptomycin Control Positive (Days)	2	4	5	9	10	10	12
Streptomycin Drug Positive (Days)	5	2	6	7	14	13	10
Result	Invalid Growth Control, Repeat Test	Invalid, Repeat Test	Resistant*	Resistant*	Susceptible	Resistant*	Invalid Growth Control, Repeat Test

\*All resistant results should be confirmed by an alternate method prior to reporting.

### LIMITATIONS

1. The VersaTREK Myco Streptomycin susceptibility test does not interpret the degree of susceptibility of the isolate being tested. It is a qualitative test.
2. The VersaTREK Myco Streptomycin susceptibility test can only be performed on the VTI or ESP instruments. Testing cannot be performed manually.
3. Use only pure cultures of *M. tuberculosis*. The user must confirm the purity of the culture, especially if the result is resistant. Cultures that contain multiple species of Mycobacteria may give erroneous results.
4. Failure to prepare the isolates per the manufacturer's recommendation may lead to inaccurate results.
5. Failure to re-hydrate the streptomycin with the proper volume of sterile water may give inaccurate results.
6. Thorough mixing of inoculated bottles is important. Failure to mix adequately may lead to false resistant results.
7. Failure to use the VersaTREK Myco GS growth supplement may lead to false susceptible results.

### QUALITY CONTROL

A new lot number or shipment of VersaTREK Myco Streptomycin should be tested with the quality control organism *M. tuberculosis* ATCC 25618 (streptomycin susceptible). Expected results are shown in Table 4. Quality control should be run once each week when susceptibility testing is performed. Observation of the proper results, within 4-13 days, indicated the reagents are ready for use in testing patient isolates.

If proper results are not observed, do not report patient's results. Repeat quality control and any patient isolates affected by the initial quality control failure. If repeat quality control does not perform as expected, do not report patient results. Do not use the product until you have contacted Technical Services.

1. Add 1 ml of VersaTREK Myco GS to VersaTREK Myco culture bottles.
2. Rehydrate lyophilized VersaTREK Myco Streptomycin with sterile water, and prepare the 5 ml to 15 ml dilution, as described under "Inoculation of VersaTREK Myco Culture Bottles for Susceptibility Testing".
3. Add 1.0 ml of the 120 µg/ml drug concentration to a VersaTREK Myco bottle. The final streptomycin concentration in the Myco culture bottle is 8 µg/ml.
4. Add 1.0 ml of the 30 µg/ml drug concentration to a second VersaTREK Myco bottle. The final streptomycin concentration in the Myco culture bottle is 2 µg/ml.
5. Add 1.0 ml of sterile distilled or deionized water to the growth control bottle.
6. Add 0.5 ml of the 1:10 dilution of culture from a positive VersaTREK Myco seed bottle to the appropriate culture bottles containing drug and to the drug-free growth control bottle.
7. After observation of expected results, use the remainder of the lot in testing specimens. If expected results are not achieved, contact Technical Support.

**Table 4: QC Strain Expected Results**

Organism	ATCC	Streptomycin
<i>M. tuberculosis</i>	25618	Susceptible

The most common causes of quality control failure are over-inoculation of the streptomycin test sets and contamination. If repeated false resistance results are obtained, the streptomycin antibiotic most likely has lost activity. The recommendation is to test with a new antibiotic stock.

**EXPECTED RESULTS**

Initial testing of 109 clinical strains of *Mycobacterium tuberculosis* was performed for susceptibility to streptomycin at 4 diverse sites using the Modified Agar Proportions (MAP) Method<sup>4</sup> and VersaTREK Myco Streptomycin Test. Both systems detected 95 (87.2%) susceptible strains. The MAP method detected 6 high resistant strains (resistant to 10 µg/ml streptomycin), and 14 low resistant strains (resistant to 2 µg/ml concentration of streptomycin). The VersaTREK Streptomycin test detected 5 high-resistant, and 11 low-resistant strains. The initial testing was performed with only one concentration (8 µg/ml) of streptomycin, which was found to be inconsistent in detecting low streptomycin resistant strains.

An additional concentration of 2 µg/ml was added to the test, and this was tested with 14 *M. tuberculosis* clinical strains, and 30 CDC challenge strains at three sites. This testing showed good correlation in detecting low resistance to streptomycin.

**PERFORMANCE CHARACTERISTICS****Table 5. Clinical Study Results with the VersaTREK/ESP Culture System II Compared to the Modified Agar Proportion Method (MAP)**

Drug	MAP			VersaTREK/ESP Culture System II				
				Susceptible Strains		Resistant Strains		
	Concentration	#S	#R	Equivalent Concentration	# Agreement	% Agreement	# Agreement	% Agreement
Streptomycin	2 µg/ml	95	14	2 µg/ml	95	100	11	79
	10 µg/ml	103	6	8 µg/ml	103	100	5	83

**Table 6. Ten Clinical (30 tests) and 30 CDC Challenge Strains (60 tests) of *M. tuberculosis* tested at two or three sites with 2 concentrations of Streptomycin (The data are represented as total number of data points tested.)**

Drug	MAP			VersaTREK/ESP Culture System II				
				Susceptible Strains		Resistant Strains		
	Concentration	#S	#R	Equivalent Concentration	# Agreement	% Agreement	# Agreement	% Agreement
Streptomycin	2 µg/ml	72	18	2 µg/ml	72	100	18	100
	10 µg/ml	77	13	8 µg/ml	77	100	13	100

**HAZARD INFORMATION**

Xn Harmful  
R22 Harmful if swallowed.  
R61 May cause harm to the unborn child.  
S37 Wear suitable gloves.  
S53 Avoid exposure – obtain special instructions before use.

**REFERENCES**

1. Stratton, C.W. 1992. The resurgence of tuberculosis. *Inf. Dis. Newsl.* 11:73-80.
2. Tenover, F.C., J.T. Crawford, R.E. Huebner, L.J. Geiter, C.R. Horsberg, Jr., and R.C. Good. 1993. The resurgence of tuberculosis: Is your laboratory ready? *J. Clin. Microbiol.* 31:767-770.
3. CDC/National Institutes of Health, Biosafety in Microbiological and Biomedical Laboratories. 4th Edition Atlanta: US Department of Health and Human Services, Public Health Service, CDC and NIH, 1999; DHHS publication no. (CDC) 93-8395.
4. Strong, B. E., and G. P. Kubica. 1981. Isolation and Identification of *Mycobacterium tuberculosis*: A Guide for the Level II Laboratory. US Publication no.81-9390. Centers for Disease Control, Atlanta, GA.

ATCC is a trademark of the American Type Culture Collection  
ESP and VersaTREK are trademarks of TREK Diagnostic Systems, Inc.

For additional information call TREK Technical Support at 1.800.642.7029 or or 1.216.351.8735.

TREK Diagnostic Systems  
982 Keynote Circle, Suite 6  
Cleveland, Ohio 44131  
USA  
800.871.8909  
[www.trekds.com](http://www.trekds.com)

TREK Diagnostic Systems  
Units 17 – 19, Birches Industrial Estate  
East Grinstead  
West Sussex, RH19 1XZ  
UK

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