INDEX, ÍNDICE, INDICE, INDEX, ÍNDICE	
2 – 8	English
9 – 15	Español
16 – 22	Italiano
23 – 29	Français
30 – 36	Português
37– 38	Symbols Símbolos Simboli Symboles Símbolos

English

VersaTREK® MYCO SUSCEPTIBILITY KIT, 7115-60

INTENDED USE

The VersaTREK Myco Susceptibility Kit is intended for qualitative *in vitro* susceptibility testing of isolated colonies of *Mycobacterium tuberculosis* with rifampin, ethambutol, and isoniazid on either the VersaTREK Instrument or the ESP Culture System II. Appropriate inoculum sources are prepared from growth on solid agar, such as Middlebrook 7H10 or Lowenstein-Jensen slants, or Middlebrook 7H9 broth.

SUMMARY AND EXPLANATION

Rapid diagnosis of tuberculosis is important for initiation of appropriate drug therapy and institution of measures to prevent the spread of this highly contagious disease. Delayed treatment, in particular of multiple-drug-resistant *M. tuberculosis* (MDR-TB) strains, can result in treatment failure and patient mortality.¹ To address these issues, clinical microbiologists are urged to use the most rapid and sensitive diagnostic tests that provide reliable results. Experts at the Centers for Disease Control (CDC) have made several recommendations regarding mycobacterial testing, including the use of both a liquid and solid medium for mycobacterial culture.²

For susceptibility testing of *Mycobacterium tuberculosis*, the VersaTREK/ESP Culture System II combines a liquid culture medium (VersaTREK Myco Broth), a growth supplement (VersaTREK Myco GS), and specific antibiotics with a detection system that automatically incubates and continuously monitors culture bottles inoculated with isolates of *Mycobacterium tuberculosis* from various specimen sources. The sponges in the VersaTREK Myco bottles provide a growth support matrix and increase the surface area exposed to headspace oxygen. The technology of the VersaTREK/ESP Culture System II is based on the detection of headspace pressure changes within a sealed bottle due to microbial growth. A special algorithm has been developed for detection of very slow growing Mycobacteria.

With the VersaTREK/ESP Culture System II, anti-tuberculosis drugs can be injected into VersaTREK Myco bottles supplemented with VersaTREK Myco GS and specific *Mycobacterium tuberculosis* isolates. An inoculum of *Mycobacterium tuberculosis*, prepared from colonies grown on solid media or from a seed bottle, is added to a drug-free control bottle, used as a positive control, and the bottles containing the specific concentrations of drug. The presence of growth is determined automatically by the instrument, but the actual determination of susceptibility or resistance is performed by a manual comparison of the drug-free positive control bottle and the drug-containing bottles. If the organism is susceptible to a drug, either no microbial growth or delayed growth will be detected. If the organism is resistant to a drug, microbial growth will be detected within a specified time and is graphically observed on the instrument by the generation of a downward consumption curve. The drugs used in this system are the primary drugs rifampin, isoniazid, and ethambutol at specific concentrations. These drugs are provided in a lyophilized form.

PRINCIPLES OF THE PROCEDURE

Susceptibility testing of *Mycobacterium tuberculosis* specimens are performed with the VersaTREK/ESP Culture System II by inoculating a suspension of organisms into VersaTREK Myco bottles that have been supplemented with VersaTREK Myco GS and with one of the three primary drugs. A VersaTREK Connector is attached to each bottle to establish a closed system for monitoring the bottle with the VersaTREK/ESP instrument. A hydrophobic membrane in the VersaTREK Connector prevents aerosolization. Clinical information is entered into the system computer and the bottles are placed in the instrument for incubation at 35 – 37 °C under stationary conditions. VersaTREK/ESP Culture System II detects mycobacterial growth by automatically monitoring (every 24 minutes) the rate of consumption of oxygen within the headspace of the culture bottle and reports that growth response with a positive signal. At the end of the specified incubation period, as determined by the drug-free control bottle for each isolate that is tested, the isolate is manually determined to be susceptible or resistant to a drug.

REAGENTS

The VersaTREK Myco Susceptibility Kit contains two bottles each of rifampin and isoniazid and three bottles of ethambutol in lyophilized form. Each bottle is to be rehydrated with 25 ml of sterile distilled or deionized water.

FORMULA UPON REHYDRATION*

(H2O)	Processed Water	25 ml
(RIF)	Rifampin	0.003% w/v
(SA)	Solubilizing Agent	1.25 % w/v
(H2O)	Processed Water	25 ml
(INH)	Isoniazid	0.0012% w/v
(SA)	Solubilizing Agent	1.25% w/v
(H2O)	Processed Water	25 ml
(EMB)	Ethambutol	0.024% w/v
(SA)	Solubilizing Agent	1.25% w/v

*Adjusted to meet performance criteria.

PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Refer to Material Safety Data Sheets for complete hazard information.
3. **WARNING! POTENTIAL INFECTIOUS TEST SPECIMEN.** Follow Universal Precautions and institutional policy in handling and disposing of infectious agents.
4. Working with *Mycobacterium tuberculosis* grown in culture, the recommendation is for Biosafety Level 3 practices, containment equipment and facilities.^{3, 4, 5}
5. The VersaTREK Connector contains a sharp recessed needle that is sterile upon first removing the seal, but is potentially contaminated upon termination of use.

Follow institutional policy for handling and disposing of blood-contaminated devices. (Do not handle in a casual manner. Sterilize prior to disposal. Dispose in an approved sharps container. Do not reuse.)

6. Visually inspect all bottles for contamination, cracks or other signs of deterioration. Do not use bottles that appear turbid or cracked.
7. Use only needle locking syringes, if using this inoculation method.
NOTE: Inoculation may be performed by opening the bottles and adding the sample using sterile pipettes.
8. Autoclave all VersaTREK Myco bottles prior to disposal.

STORAGE

VersaTREK Myco Susceptibility Kit: Store at 2-8°C. Protect from light. Upon rehydration, store at 2-8°C and use within 4 weeks, or freeze at -20°C and use within 3 months. Once thawed, do not refreeze.

PROCEDURE: SUSCEPTIBILITY TESTING OF *MYCOBACTERIUM TUBERCULOSIS* USING THE VersaTREK/ESP CULTURE SYSTEM II

Materials Provided:

VersaTREK Myco Culture Bottles	7111-42
VersaTREK Myco GS	7112-42
VersaTREK Myco Susceptibility Kit	7115-60
VersaTREK Connectors	7150-44

Materials Required but Not Provided:

Sterile distilled or deionized water
Middlebrook 7H10 or other mycobacterial agar or egg-base medium
Autoclave
CO₂ incubator, 37°C
Sterile loops
Sterile tuberculin syringes with permanently attached needles
Sterile Pipettes or Pipettor with Sterile Pipette Tips (optional)
Mycobactericidal disinfectant (e.g. AmphyI™, OMNI™, etc.)
Alcohol swabs
Vortex mixer
Biological safety cabinet
Microscope
Materials for staining slides
Quality Control organisms: *M. tuberculosis* ATCC® 25618, *M. tuberculosis* ATCC 35838, *M. tuberculosis* ATCC 35822, and *M. tuberculosis* ATCC 35837

Inoculation of VersaTREK Myco Culture Bottles for Myco Susceptibility Testing

1. Rehydrate the lyophilized drugs by aseptically adding 25 ml sterile distilled or deionized water to each drug vial. Additional drug concentrations should be made for the isoniazid and ethambutol vials as instructed in Table 1.

Table 1: Anti-Tuberculosis Drug Preparation

Drug and Final Concentration in Myco Bottles	Amount of Diluent Added to Lyophilized Bottles	Concentration after Rehydration	Additional Dilutions	Concentration to be Injected into Test Bottles	Number of Tests/Vial	Number of Tests/Kit
Rifampin 1 µg/ml	25 ml	30 µg/ml	No	30 µg/ml	50	100
Isoniazid 0.4 µg/ml Isoniazid 0.1 µg/ml	25 ml	12µg/ml	No 5 ml diluted drug plus 15 ml water	12 µg/ml 3 µg/ml	40 40	80 80
Ethambutol 8 µg/ml Ethambutol 5 µg/ml	25 ml	240µg/ml	No 10 ml diluted drug plus 6 ml water	240 µg/ml 150 µg/ml	30 32	90 96

NOTE: There is sufficient volume to perform 80 tests per kit. The above dilutions provide an excess of some of the drugs in some concentrations. Discard excess amounts and do not mix components of this kit with different lot numbers of other kits.

2. The inoculum of a specimen for susceptibility testing can originate from two sources. Prior to setting up a susceptibility test on any isolate, a determination should be made that the isolate is *Mycobacterium tuberculosis* and that it is a pure state (contaminant-free):

a. McFarland equivalent cell suspension - Isolates growing on an appropriate mycobacterial agar medium (Middlebrook 7H10, Middlebrook 7H11, Lowenstein-Jensen) are scraped using a sterile loop and placed in a tube containing approximately 9 to 10 ml sterile distilled or deionized water and glass beads. The tube is capped, vortexed, and allowed to settle for 15 minutes. Approximately 5 ml of the top layer is removed, placed in a sterile test tube, and adjusted to a No. 1 McFarland equivalent (approximately 3.0×10^7 CFU/ml) using sterile distilled or deionized water. A 1:10 dilution of the adjusted sample is prepared using sterile distilled or deionized water. For the final inoculum, the organism is diluted 1:10.

b. VersaTREK Myco Seed Bottle - Prepare a No. 1.0 McFarland equivalent created from an agar or liquid source (growth in Middlebrook 7H9) of the test isolate. Aseptically inject 0.5 ml of this dilution into a VersaTREK Myco culture bottle supplemented with 1.0 ml VersaTREK Myco GS. A sterile VersaTREK Connector is placed on the bottle and the bottle is incubated in the instrument until a positive signal is observed (approx. 10^6 CFU/ml). Once a positive signal is reached, the VersaTREK Connector is removed and the bottle is vortexed for 1-2 minutes. The seed bottle inoculum must be used within 72 hours of the positive instrument signal. Dilute the seed bottle 1:10.

3. Label each VersaTREK Myco culture bottle with Patient Identification or QC identifier. For Myco Susceptibility testing, 6 bottles will be used per isolate.

4. Disinfect the culture bottle stoppers with alcohol and aseptically add the following components into each culture bottle by injecting through the septum with a needle and locking syringe as described in Table 2. NOTE: Alternatively, these additions may be performed by opening the bottles and adding the reagents and sample using sterile pipettes. Care must be taken to maintain aseptic technique. The stopper and cap must be resealed properly to guarantee a sealed system. If the stopper and cap are not seated properly, a leak bottle error will occur.

Table 2: Reagent and Inoculum Additions to the Test Procedure

Label Myco Bottle with final concentration and patient identification number	Volume of diluted GS to add	Volume of diluted drug prepared in Table 1 to add	1:10 dilution from either source of prepared inoculum
1 µg/ml Rifampin	1 ml	0.5 ml Rifampin at 30µg/ml	0.5 ml
0.4 µg/ml Isoniazid	1 ml	0.5 ml Isoniazid at 12 µg/ml	0.5 ml
0.1 µg/ml Isoniazid	1 ml	0.5 ml Isoniazid at 3 µg/ml	0.5 ml
8 µg/ml Ethambutol	1 ml	0.5 ml Ethambutol at 240 µg/ml	0.5 ml
5 µg/ml Ethambutol	1 ml	0.5 ml Ethambutol at 150 µg/ml	0.5 ml
Control	1 ml	0.5 ml sterile distilled or deionized water	0.5 ml

5. Perform a purity check on the *Mycobacterium tuberculosis* inoculum by subculturing the organism suspension on both a Middlebrook 7H10 or 7H11 plate and a solid medium that will support the growth of most facultative bacteria, such as Blood Agar (BAP). Incubate the BAP for at least 18 hours in a conventional incubator and the Middlebrook agar for 4 weeks at 37°C.
 6. Wipe the culture bottles and stoppers with mycobactericidal disinfectant.
 7. Mix all bottles by inverting several times.
 8. Remove the seal from the bottom of a VersaTREK Connector. Press the needle end of the VersaTREK Connector through the culture bottle stopper. NOTE: Do not invert the culture bottle while the VersaTREK Connector is on the bottle. Fluid in the needle may interfere with the pressure readings for the bottle.
 9. Label each VersaTREK Myco bottle with Patient ID, Drug ID, and concentration. Recommended format for the accession number is "Access" XY, where "Access" equals the accession number for a particular isolate's susceptibility test, X equals the drug identifier, and Y equals the drug concentration. This labeling strategy will allow the user to retrieve the test results for each isolate group of bottles. You can use the "Access" portion of the number to search for and display the test results. Refer to the VersaTREK/ESP Culture System II Operator's Manual for details on bottle entry.
- | | | |
|----------|---------------------|-----------------------------|
| Example: | Control | Accession number = 1234C |
| | Rifampin 1.0 µg/ml | Accession number = 1234R1 |
| | Isoniazid 0.4 µg/ml | Accession number = 1234I0.4 |
10. Place the culture bottle with VersaTREK Connector into the instrument.
 11. Detection of growth is determined by the instrument and the time to detection can be obtained by reviewing the graph of the test result. The test should be completed in 13 days or less based on the drug-free control bottle.

RESULTS

For susceptibility testing using the VersaTREK/ESP Culture System II, an isolate is interpreted as being susceptible or resistant to a drug based on the following (round all times of detection to the nearest whole day):

1. After the drug-free bottle has become positive, record the time to detection (this is found on the graph of the test bottle).
2. Three additional days of incubation after the drug-free control has signaled positive, search and record the results from the other test bottles for the isolate.
3. **Susceptible:** An isolate in a drug-containing bottle that has not signaled within the three-day detection period of the drug-free control is determined to be susceptible at that drug concentration. Growth after the three days, rounded to the nearest whole number, of the growth control bottle signaling positive should be disregarded because it is not considered to be an indication of resistance. To round to the nearest whole number, if time \leq XX.49 is round down and if time \geq XX.50 is round up.
4. **Resistant:** Time to detection of drug bottle is \pm three days (rounded to the nearest whole number) of the time to detection of the drug-free control bottle (rounded to the nearest whole number). To round to the nearest whole number, if time \leq XX.49 is round down and if time \geq XX.50 is round up.

NOTE:

- No bottle should flag positive before the third day. Bottles that flag before three days are either contaminated or the inoculum is too high. The test should be repeated.
- A test bottle may signal before the control bottle – see Table 4 example 2.

Bottles that signal positive should be stained for AFB (Acid Fast Bacilli) and the plates from the initial inoculum preparation that was subcultured to a Middlebrook 7H10 or 7H11 agar plate should be observed for growth at this time and for up to 4 weeks. If non-AFB smear organisms appear to be present, a subculture to a non-selective plate (e.g. BAP) should be done and incubated for at least 18 hours in a conventional incubator at 35°C to verify purity of the isolate. If it is necessary to repeat the testing, the subculture from the original inoculum for purity verification may be used. Test observations that determine the validity of the test are shown in Table 3. An example of results for determination of susceptibility is shown in Table 4.

TABLE 3: Test Observations and Test Assessment

Test bottle	Result obtained	Possible cause	Action
Growth control	<3 days from inoculation	Possible contaminant or too heavy inoculum; check purity plates and graphs.	REPEAT TEST
Growth control*	>10 days from inoculation	Non viable organism, too light inoculum	REPEAT TEST
Growth control	Growth detected between 3 and 10 days after inoculation		Valid test - evaluate drug bottles
Purity check	No contaminants observed by at least 18 h on Blood Agar plate		Valid test
	Non-AFB organisms observed	Culture was not pure	REPEAT TEST
	No colonial variation observed on Middlebrook agar plates		Valid test- No action
	Colonial variation observed on Middlebrook agar plates	More than one species of mycobacterium is present	REPEAT TEST to verify results; notify physician that final results are pending

*The upper range of 10 days detection time for the growth control was established by clinical and in-house data.

Be sure to check the subculture results on the mycobacterial agar to verify that only one species of mycobacterium was present, especially in the case of a resistant result.

Always review the isolate's graph before reporting. Refer to the Operator's Manual for graph analysis. If graphs for tests that have signaled positive look abnormal, contact TREK Technical Support at 1.800.642.7029 or 1.216.351.8735.

Table 4: Examples of Results for Calculation of Susceptibility of *M. tuberculosis* Strains to VersaTREK Myco Susceptibility Kit Antibiotics

Example 1

Drug/ Conc. (µg/ml)	Rifampin	Isoniazid		Ethambutol	
	1	0.1	0.4	5	8
Control Positive	2				
Day Positive	4	5	7	5	No Growth
Result	Invalid Growth Control, Repeat Test				

Example 2

Drug/ Conc. (µg/ml)	Rifampin	Isoniazid		Ethambutol	
	1	0.1	0.4	5	8
Control Positive	5				
Day Positive	3	7	9	9	No Growth
Result	Resistant	Resistant	Susceptible	Susceptible	Susceptible

Example 3

Drug/ Conc. (µg/ml)	Rifampin	Isoniazid		Ethambutol	
	1	0.1	0.4	5	8
Control Positive	7				
Day Positive	8	11	No Growth	10	No Growth
Result	Resistant	Resistant	Susceptible	Resistant	Susceptible

Example 4

Drug/ Conc. (µg/ml)	Rifampin		Isoniazid		Ethambutol	
	1	0.1	0.4	5	8	
Control Positive	11					
Day Positive	9	11	No Growth	10	No Growth	
Result	Invalid Growth Control, Repeat Test					

LIMITATIONS

1. The ability of the VersaTREK/ESP Culture System II to detect resistance to low concentrations of Ethambutol with *Mycobacterium tuberculosis* is unknown because testing at the concentration of 5 µg/ml was not done during the test evaluation at the clinical sites.
2. The VersaTREK Myco Susceptibility Kit should not be used with isolates obtained directly from a positive specimen bottle (Primary Bottle).

EXPECTED RESULTS

109 Clinical strains of *Mycobacterium tuberculosis* were tested for susceptibility to Rifampin (RIF), Ethambutol (ETH), and Isoniazid (INH) at 4 diverse sites using the Modified Agar Proportions (MAP) Method⁷ and the VersaTREK Myco Susceptibility Test System. Both methods detected 87 (79.8%) pan-susceptible strains. The MAP method detected 11 mono-resistant INH strains (10.1%), 9 Multiple Drug Resistant (MDR) Strains (8.2%), and 2 strains (1.8%) that were resistant to all three drugs. One strain was resistant to INH and ETH. The VersaTREK/ESP Culture System II detected 11 INH (10.1%) and 1 RIF (0.9%) monoresistant strains, 8 MDR strains (7.3%), one strain resistant to INH and ETH, and two strains (1.8%) that were resistant to all three drugs tested. According to the WHO/IUALTD Global Project on Anti-Tuberculosis Drug Resistance Surveillance⁶, the most common resistance patterns in the United States are mono-resistance to Isoniazid (INH) and Streptomycin (STR), followed by Multiple Drug Resistance (MDR). In the surveillance report, 87.1% of the strains tested were susceptible. The percent resistance to multiple drugs was 4.5%, of which 2% were MDR. These data are in line with the study results reported with VersaTREK.

QUALITY CONTROL

Each new VersaTREK Myco Susceptibility Kit should be tested with the quality control isolates listed in Table 5 prior to patient testing. Quality control strains should be tested the same as clinical isolates in the procedure for "Inoculation of VersaTREK Myco Culture Bottles for Myco Susceptibility Testing." Results should be read as recommended in the Results section of this package insert. The following chart provides the expected results for the recommended quality control strains. If the expected results are not achieved, do not perform or report patient results using these lots. Contact TREK Technical Support at 1.800.642.7029 or 1.216-351.8735.

Table 5: QC Strains and Expected Results

Organism	ATCC	Rifampin	Isoniazid	Ethambutol
<i>M. tuberculosis</i>	25618	Susceptible	Susceptible	Susceptible
<i>M. tuberculosis</i>	35838	Resistant	–	–
<i>M. tuberculosis</i>	35822	–	Resistant	–
<i>M. tuberculosis</i>	35837	–	–	Resistant

PERFORMANCE CHARACTERISTICS

A comparison between VersaTREK/ESP Culture System II and the Agar Proportions Method⁷ was performed with a total of 109 clinical strains of *Mycobacterium tuberculosis* from four test sites. The test sites included an urban hospital with large HIV patient populations and high Mycobacteria test volume, two national TB referral hospitals with high volume of referral cultures, and a large university hospital. Testing was performed using three of the primary drugs for *M. tuberculosis* susceptibility testing: Rifampin, Isoniazid, and Ethambutol. Both VersaTREK Myco broth seeded and the McFarland standard methods of inoculation preparation were evaluated with similar results; therefore only one set of data is presented. Data from two strains were eliminated because of contamination; one strain from the VersaTREK/ESP Culture System II and one from the agar proportions method. The testing was performed on one concentration of Rifampin and two concentrations of Isoniazid and Ethambutol. Table 6 presents the agreement and 95% confidence intervals (CI) of the VersaTREK/ESP Culture System II with the agar proportions test results.

**Table 6: Clinical Study Results with the VersaTREK/ESP Culture System II
Compared to the Agar Proportions Method**

Drug	Agar Proportion Method			VersaTREK/ESP Culture System II				
	Concentration	# S	# R	Equivalent Concentration	Susceptible Strains		Resistant Strains	
					#Agreement	%Agreement (95% CI)	#Agreement	%Agreement (95% CI)
Rifampin	1 µg/ml	99	10	1 µg/ml	99	100 (96-100)	10	100 (69-100)
Isoniazid	0.2 µg/ml	85	24	0.1 µg/ml	85	100 (96-100)	23	95.8 (52-100)
	1 µg/ml	94	15	0.4 µg/ml	94	100 (96-100)	15	100 (79-99)
Ethambutol	5 µg/ml	104	4	5 µg/ml*	NT		NT*	
	10 µg/ml	105	3	8 µg/ml	105	100 (96-100)	3	100 (29-100)

* Repeat testing using the VersaTREK/ESP Culture System II for Ethambutol at 5 µg/ml gave a resistant result for the one isolate that was susceptible at 10 µg/ml but resistant at 5 µg/ml when tested with the Agar Proportions Method.
NT = Not Tested

HAZARD INFORMATION

Rifampin

Xn Harmful
R22 Harmful if swallowed.
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S36/39 Wear suitable protective clothing and eye/face protection.

Isoniazid

Xn Harmful
R22 Harmful if swallowed.
S20 When using do not eat or drink.
S36 Wear suitable protective clothing.
S46 If swallowed, seek medical advice immediately and show this container or label.

Ethambutol

Xn Harmful
R22 Harmful if swallowed.
S20 When using do not eat or drink.
S24/25 Avoid contact with skin and eyes.
S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

REFERENCES

1. Stratton, C.W. 1992. The Resurgence of Tuberculosis. *Int. Dis. Newsl.* 11:73-80.
2. Tenover, F.C., J.T. Crawford, R.E. Huebner, L.J. Geiter, C.R. Horsburg, Jr., and R.C. Good, 1993. The Resurgence of Tuberculosis: Is Your Laboratory Ready? *J. Clin. Microbiol.* 31:767-770.
3. CDC Guidelines for Preventing the Transmission of Mycobacterium tuberculosis in Health-care Facilities. 1994. *MMWR* 1994, 43 (No. RR-13)
4. CDC/National Institutes of Health, Biosafety in Microbiological and Biomedical Laboratories. 4th Edition Atlanta: US Department of Health and Human Services, Public Health Service, CDC and NIH, 1999; DHHS publication no. (CDC) 93-8395.
5. Gilchrist, MJR, Biosafety Precautions for Airborne Pathogens. In: Flemng, D.O., J.H. Richardson, J.J. Rylis, D. Vesely, eds., *Laboratory Safety: Principles and Practices*, 2nd Edition. Washington, D.C. American Society for Microbiology Press, 1991; 4:207-241.
6. Anti-Tuberculosis Drug Resistance in the World. The WHO/IUALTD Global Project on Anti-Tuberculosis Drug Resistance Surveillance. 1994-1997.
7. Strong, B.E., and G.P. Kubica. 1981 Isolation and identification of Mycobacterium tuberculosis: A Guide for the Level II Laboratory. US Publication no. 81-9390. Centers for Disease Control, Atlanta, GA.

For additional information, contact TREK Technical Support at 1.800.642.7029. Outside the U.S., call 1.216.351.8735.

TREK Diagnostic Systems
982 Keynote Circle, Suite 6
Cleveland, Ohio 44131
USA
800.871.8909
www.trekds.com

TREK Diagnostic Systems
Units 17 – 19, Birches Industrial Estate
East Grinstead
West Sussex, RH19 1XZ
UK
44.1342.318777

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